TUNTURI®

USER MANUAL

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SC30 SMART SCALE

14TUSSP002

TUNTURI®	
SC30 Smart Scale	
SKU: 14TUSSP002	
INPUT: 3V DC (2x AAA Battery)	
BT Address: Tunturi SC30	
MAX USER WEIGHT: 180 KGS	
HOME USE	
Serial nr:	
CE RoHS FC Made in China	

FEEL BETTER **EVERY DAY**

TABLE OF CONTENT

English	3
Product specifications	4
Product features	4
14 body measurements	5
Caution	6
Troubleshooting	7
How to weigh	8
Operating instructions	12
Maintenance	12
Declaration of the manufacturer	13
Disclaimer	13
Nederlands	14
Productspecificaties	15
Productkenmerken	15
14 lichaamsmetingen	16
Waarschuwing	18
Probleemoplossing	18
Hoe weeg je	19
Gebruiksaanwijzing	19
Onderhoud	23
Verklaring van de fabrikant	24
Disclaimer	24
Deutsch	25
Produktspezifikationen	26
Produktmerkmale	26
14 Körpermaße	27
Vorsicht	28
Fehlerbehebung	29
Wie man wiegt	30
Bedienungsanleitung	30
Wartung	34
Herstellererklärung	35
Disclaimer	35
Français	36
Spécifications du produit	37
Caractéristiques du produit	37
14 Mesures corporelles	38
Mise en garde	39
Dépannage	40
Comment peser	41
Instructions d'utilisation	42
Maintenance	45
Déclaration du fabricant	46
Disclaimer	46

Thank you for choosing the Tunturi Smart Scale! This scale is your personal health assistant. It uses bio-electrical impedance analysis (BIA) technology to offer you the data you need to keep track of your personal health metrics: BMI (Body Mass Index), Body Fat %, Body Water, Muscle Mass, Bone Mass, Protein, and much more!

Note:

The maximum user weight is 180 kg.

PRODUCT SPECIFICATIONS


EN

Dimensions (LWH)	26 x 26 x 1.8 cm 10.2 x 10.2 x 0.7 inch
Weight	1 kg 2.2 lb
Batteries	2x 1.5V AAA
Max. user weight	180 kg 396 lb
Division	0.1 kg 0.22 lb
Unit	kg lb

ACCURACY RANGE

50 kg ±0.6%	100 kg ±0.4%	180 kg ±0.3%
110 lb ±0.6%	220 lb ±0.4%	396 lb ±0.3%

PRODUCT FEATURES

1. High Sensitivity Electrodes
2. Durable Tempered Glass
3. LCD display

4. Anti-slip pad
5. Battery compartment
6. Switch button: kg/lb
7. 2x 1.5V AAA Batteries

14 BODY MEASUREMENTS

EN

WEIGHT

BMI

Body Mass Index

BODY FAT

MUSCLE RATE

FAT-FREE BODY WEIGHT

SUBCUTANEOUS FAT

VISCERAL FAT

BODY WATER


SKELETAL MUSCLE

MUSCLE MASS

BODY AGE

PROTEIN

BMR

Basal Metabolic Rate

BONE MASS

CAUTION**EN****Safety precautions:**

When using digital scale, certain basic safety precautions must be followed, including the following:

1. The Bluetooth scale should not be used to diagnose or treat any medical condition. You should consult a doctor when undertaking any diet or exercise program.
2. If you are pregnant or have a pacemaker or other internal device, you should not rely on the data provided by this product.
3. All data is for reference and is not a substitute for advice from a licensed health care professional.
4. Check the device before each use. Do not use the device if damaged. The continuous use of a damaged unit may cause injury or improper results.

Points of attention:

1. Scale is a high-precision measuring device. Never jump or stamp on the scale or disassemble. Please handle the scale carefully to avoid breakage.
2. Place the scale on a hard, flat surface during use. Using the scale on a soft surface like carpet or linoleum may affect the accuracy of the scale readings.
3. To ensure weighing accuracy, please gently step on the scale to wake it up and wait at least 10 seconds before weighing yourself on the scale. You must also wake up the scale if it has been moved.
Take your measurements at the same time each day to ensure the most accurate results.
4. When not in use, please remove the batteries and place the device in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the equipment.
5. Before using the scale, make sure you have entered all your personal data through the Fitdays app.
6. Use the scale while barefoot. The scale cannot measure body fat percentage if you're wearing shoes or socks.
7. Make sure your feet are dry before weighing yourself.
8. This scale is not waterproof; never immerse the scale in water. To clean the surface, use a damp cloth or glass cleaner. Do not use soap or other chemicals.

CAUTION

EN

Battery safety:

1. If the scale is not used for long periods, it's recommended to remove the battery to avoid damage to the scale due to possible battery leakage.
2. Keep the battery out of reach of children.
3. Do not mix old and new batteries, with different compositions or of different brands to prevent possible leakage or explosion.
4. Do not heat or deform the batteries or expose to fire.
5. Waste batteries should not be disposed of with household waste.
6. Please check with your local authority for battery recycling advice.

TROUBLESHOOTING

888.8

FUNCTIONING NORMALLY

No display

CHANGE BATTERY

[

CALIBRATING

Lo

CHANGE BATTERY

Err



OVERLOAD

If the scale malfunctions, first check the batteries and replace if needed.



HOW TO WEIGH

EN

For body weight: Place the scale on an even surface and wait at least 10 seconds before stepping on the scale. The reading should display 0.0 kg before stepping on.

For body fat and other body composition analysis: Stand straight on the scale and uniformly touch the electrodes with bare and dry feet, as shown below.

Conductive areas

If you want to switch between kg and lb, press the button on the back of the product. You can see the kg/lb indication changing on the LCD screen.

OPERATING INSTRUCTIONS

The Bluetooth Scale can connect to your mobile phone (Android 6.0+ & IOS 8.0+) or iPad via Bluetooth 4.0. Download the "Fitdays" app free from the App Store or Google Play.

Please note that the instructions and screenshots can change with app updates.


OPERATING INSTRUCTIONS

EN


DOWNLOAD APP AND INSTALL BATTERIES

Download "Fitdays" via App Store or Google Play.


Register your ID, then add personal data and confirm "ok".


Install the original battery set. Please note the batteries offered are trial versions. We recommend replacing them within three months and removing the batteries when not used for a long time.

PAIRING THE SCALE WITH YOUR SMARTPHONE


Turn on Bluetooth on your smartphone. Location also needs to be turned on for Android 6.0 or higher.

Press [Device] and Device Management to enter measuring page and find device.

Power on the scale, press[+] to select "Tunturi SC30" to add the device.


Come back to measure page. Step on the scale and keep standing for 5-10 seconds with bare feet. The weighing has finished after the analysing data shows.


OPERATING INSTRUCTIONS

EN


PAIRING THE SCALE WITH YOUR APPLE WATCH

Open the Fitdays app on your Apple Watch.



If it says to open the app on your phone you have to do this step.

You can now analyze your weight and subcutaneous fat.

VIEW REPORTS, SET REFERENCE DATA, SHARE PROGRESS

Click "Charts" to view your progress history. Enter it into the calendar to form a list and select a dataset or all data you want to delete. You can also select different data to contrast and share with others.

May 2023							+
Sun	Mon	Tue	Wed	Thu	Fri	Sat	
30	01	02	03	04	05	06	
07	08	09	10	11	12	13	
14	15	16	17	18	19	20	
21	22	23	24	25	26	27	
28	29	30	31	01	02	03	
04	05	06	07	08	09	10	
16:37							
81.0 kg		23.9		22.1 %			
Weight		BMI		Body Fat			
15:32							
80.9 kg		23.9		22.1 %			
Weight		BMI		Body Fat			




OPERATING INSTRUCTIONS

EN

ADD/DELETE USERS & BASIC SETTINGS

Under "Account" you can add a new user by selecting "Add User". To delete an account, swipe the account to the left and select "Delete".

Under "Account" you can go to the "Settings" page. Within this menu you can customize the app to your preference.

SYNCING FITNESS APPS

**Apple Health**

1. Go to Apple Health, Account top right, Under Privacy, choose apps.
2. Select FitDays, Turn on all categories to allow the Fitdays app to work.

**Google Fit**

1. Within the FitDays app, Select Account and go to Settings
2. Select Google Fit and Turn on Google Fit.

OPERATING INSTRUCTIONS

EN


1. Enter the "Settings" page, and select Fitbit.
2. Turn on the Fitbit button, and login to your Fitbit account.
3. Select "Allow" to access Fitbit App.

MAINTENANCE

If the scale has been moved or flipped up side down, It must be recalibrated to ensure accurate results.

1. Place the scale on a hard, flat surface.
2. Step on the scale with one foot until the digits appear on the display, then, step off.
3. The scale will show "0.00", indicating that the calibration process is complete.

DECLARATION OF THE MANUFACTURER

EN

Tunturi New Fitness BV declares that the product is in conformity with the following standards and directives:

RoHS Directive (2011/65/EU)

Battery Regulation (EU) 2023/1542

Hereby, Tunturi New Fitness BV declares that the radio equipment type smart scale is in compliance with Directive 2014/53/EU.

The full text of the EU declaration of conformity is available at the following internet address: <https://doc.tunturi.com>

Tunturi New Fitness BV
Purmerweg 1
1311 XE Almere
The Netherlands

DISCLAIMER

© 2023 Tunturi New Fitness BV

All rights reserved. The product and the manual are subject to change. Specifications can be changed without further notice.

Tunturi New Fitness B.V. / P.O. Box 60001 / 1320 AA Almere / The Netherlands / www.tunturi.com

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